Chela Liquid

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EN Chemical reaming and cleaning of root canals

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✉ info@famedical.net

🌐 www.famedical.net

For chemical reaming and cleaning of root canals and exposure of root canals orifices. Chela Liquid is an effective EDTA solution for chelating, cleaning and debriding of root canals. EDTA reacts with calcium in canal walls helping in reaming process during mechanical instrumentation.

Chela Liquid contains benzalconium chloride providing the antibacterial effect.

•• MAIN ADVANTAGES

- Chemical expansion of root canals due to the dissolution of calcium components of the canal walls and dentin sawdust during mechanical processing.
- The material is based on EDTA: sodium salt (concentration – 15%).
- Quick dissolution of calcium of dentin plugs.
- Penetration even into the thinnest canal.
- Contains benzalkonium chloride which provides antibacterial effect of the material.
- The material minimizes the risk of tool breakage inside the canal.
- Facilitates the removal of broken tools, silver and gutta-percha pins, etc. from the canal.

•• COMPOSITION

- EDTA 15 g / Benzalconium chloride 0.3 g / Excipient Q.S. 100 g.

•• INDICATIONS

- Chemical expansion of narrow and hard-to-reach root canals with simultaneous disinfection.
- Detection of canal orifice.
- Unsealing of previously sealed canals.

•• DIRECTIONS FOR USE

- Introduce Chela Liquid into the root canal using a metal cannula or cotton pellet and leave it for 1 – 2 minutes. Low viscosity solution penetrates into the canal readily with fire. Remove softening dentin using appropriate endodontic instrument.
- Repeat the chemical and mechanical reaming operations till suitable canal preparations is obtained.
- In especially difficult cases leave the liquid in the cavity on a cotton pellet for 1 – 2 days capped with a provisional filling.

•• STORAGE AND SHELF LIFE

- The Chela Liquid bottle should be tightly closed immediately after use. Protect from direct light.
- Shelf life at room temperature is 2 years.
- Refrigeration will further increase shelf life. DO NOT FREEZE!

•• WARRANTY

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty. The customer is responsible for determination of suitability of this product for user's application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

•• LIMITATION OF LIABILITY

The manufacturer's liability is limited by only cases stipulated by direct legislation of the country.

•• RECYCLING

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

•• PACKING

- Chela Liquid – bottle 9 ml.

For professional use in dentistry only!

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